



Analysis of Sialic Acids by HPAEC-PAD

MATERIALS

- 4N Acetic Acid (HOAc)
- 50% NaOH (Fisher, SS254-1)
- 50% Isopropanol
- Sodium Acetate Trihydrate (NaOAc) (VWR, JT3462-1)
- Heat Block
- Glass Culture Tube, 13 x 100 mm, with Teflon lined Cap (Pyrex, 9826-13)

Preparation of 1L of 100 mM NaOH with 5 mM NaOAc:

- 0.7 g of Sodium Acetate Trihydrate
- 5.2 mL of 50% NaOH

Dissolve the Sodium Acetate Trihydrate in Milli-Q water that has been filtered through a 0.2 μ m filter. Then add the 50% NaOH and bring the volume to 1 L. Thoroughly mix the solution and sparge the solution with Helium for 15 min to remove dissolved gases.

Preparation of 1L of 100 mM NaOH with 250 mM NaOAc:

- 34 g of Sodium Acetate Trihydrate
- 5.2 mL of 50% NaOH

Dissolve the Sodium Acetate Trihydrate in Milli-Q water that has been filtered through a 0.2 μ m filter. Then add the 50% NaOH and bring the volume to 1 L. Thoroughly mix the solution and sparge the solution with Helium for 15 min to remove dissolved gases.

Preparation of 4 mL of 4 N Acetic Acid:

- 920 μ L of Glacial Acetic Acid

In a glass vial dissolve 920 μ L of glacial acetic acid in 3080 μ L of Milli-Q water. Store at 4 $^{\circ}$ C

PROCEDURE:

1. Dissolve 500 μ g of sample and dissolve in 200 μ L of Milli-Q water and place in a glass reaction tube. If the analysis is to be performed on a liquid sample use no more than 200 μ L for analysis.
2. To the sample add an equal volume of 4N HOAc, the final concentration of HOAc should be 2N, and tightly cap the reaction tube.

NOTE: Avoid salt buffers such as Tris, NaCl, and NaOH as this can alter the separation and retention times of the sample components

3. Hydrolyze the sample on a heating block at 80 $^{\circ}$ C for 3hr.



4. Once hydrolysis is complete remove the sample tubes from the heating block and allow them to cool to room temperature. Centrifuge the samples at 2000 rpm for 2 min to bring down any condensation on the sides and cap of the tube.
5. Next, evaporate the samples under a flow of dry nitrogen while applying a low amount of heat to the sample tubes.
6. When the samples are dry add 100 μ L of 50% IPA, vortex and evaporate the samples a second time.
7. Resuspend the dried samples in Milli-Q water. Use a volume appropriate for the amount of sample that will be injected into the HPLC.

NOTE: If this is an unknown sample resuspend in 100 μ L of Milli-Q water and inject up to 50% of the sample. Never inject 100% of the sample unless you know that it is a very small amount of material as a large amount of sample will hamper accurate quantitation of monosaccharides.

NOTE: It is essential to remove any protein from the sample using spin filtration to prevent contamination and blockage of the column, guard column and their inlet frits with protein.

HPAEC-PAD ANALYSIS OF SIALIC ACIDS:

Using Dionex ICS-3000

- Colum:
 - Dionex CarboPac PA1 column 4 mm x 250 mm, 4 μ m, with a 4 mm x 50 mm Guard
- Solvents:
 - A: Water
 - B: 100 mM NaOH with 5 mM NaOAc
 - C: 100 mM NaOH with 250 mM NaOAc
 - Initial conditions of 16% B at 1.0 mL/min
- Pulsed Amperometric Detector:
 - Waveform: Standard Quad
- Gradient Settings:

Time	%A	%B	%C
0	80	10	10
20	30	10	60
25	10	10	80
30	10	10	80
31	80	10	10
50	80	10	10