



## **Permethylation of N- and O-Glycans**

### **MATERIALS**

- Mortar and Pestle
- Dimethyl Sulfoxide (Sigma, 276855)
- NaOH Pellets
- Iodomethane (Sigma, 289566)
- Chloroform (Sigma, 288306)

### **PROCEDURE:**

NOTE: In order to perform permethylation the sample needs to be free of water. It is best to place the sample in a glass reaction tube and lyophilize it to thoroughly dry the sample.

1. Dissolve the dried glycans sample in 500  $\mu$ L of DMSO, add a stirbar, and stir for 15-30 min, while the sample tube is capped.
2. Place 2 pellets of NaOH per sample (10-20  $\mu$ g) in a mortar and crush it into a powder using the pestal. Add DMSO to the freshly ground NaOH, making a slurry.
3. While stirring with the pipette tip take up 500  $\mu$ L of the slurry and add it to the sample tube and recap.
4. Using a gas-tight syringe add 200  $\mu$ L of methyl iodide to the sample and recap. Allow the sample to stir for 30 min. The reaction is complete when the mixture is milky white in appearance. Add 50  $\mu$ L of methyl iodide to the reaction mixture and allow it to react for another 15-20 mins. (This addition should be performed in a chemical hood.)

NOTE: Methyl iodide should be kept in a container with desiccant at 4 °C. Monitor the temperature of the reaction and place the sample tubes on ice for 1 min when the sample becomes warm and then return it to stirring.

5. Place the sample tubes on ice and add 1 mL of cold Milli-Q water to stop the reaction. The contents of the tube will bubble as water is added.
6. Extract the permethylated glycans by adding 1 mL of chloroform to the sample and vortexing. Centrifuge the sample at 2000 rpm for 2 min.



7. After centrifugation remove the water layer and discard. Again add water and vortex the sample and then allow it to sit for 2 min for the phases to clearly separate. Remove the water layer and discard. Repeat this step 3-4 times.
8. Remove the chloroform using a stream of dry nitrogen.
9. Once dried the sample is ready for MALDI-TOF analysis or conversion to PMAA for linkage analysis.